




Coil Extends to 7'  
Useful Length



General: Meets or exceeds the requirements listed in ASTM D 6867-03  
 Temperature: -30 °F to 300 °F to retain configuration  
 Flammability: UL 94 VO rated  
 FDA: Material approved for use in food contact applications in compliance with FDA regulation: 21 CFR 177.1550  
 USP CLASS VI: Material meets the requirements of USP Class VI  
 NSF: Material certified to NSF 51, NSF 61 and NSF 372

Designer			TITLE: PFA Retractable Coil, 2" Coil I.D. X 2 3/4" O.D. Tubing is 3/8" O.D X 1/4" I.D., (.0625 Wall Thickness)
Review			
Engineer			
 <b>MARQUEST SCIENTIFIC</b> <small>Fluid Handling Products</small>			DWG. NO. <b>C3-PFA</b>
UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES			SCALE: N/A    Format: SD    SHEET 1 OF 1