



Notes:

*For this model, joints A, B, C & D are permanent, non-rotating, either solvent or fusion welded.

**Contact factory for options including threaded joints, swivel, and removable barb tip.

Item	Part Number	Description	Qty.
1	TPWKL-375	Top Works Sub-Assembly	1
2	NV-S-PTFE	Valve Stem Seal	1
3	LG-CWS-BUNA	Casual Water Seal	
4	LR-375	Nut Ring, Mounting	1

DRAWN	NAME	DATE
CHECKED	P.S.	
ENG APPR.	C.W.	
MFG APPR.		
Q.A.		
COMMENTS:		

Marquest Scientific, Inc.
2950 Airway Avenue, Bldg. C5, Costa Mesa, CA 92626

**Repair Parts List
for Lab Faucet LG-DN-R**

PROPRIETARY AND CONFIDENTIAL
THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF MARQUEST SCIENTIFIC, INC.. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF MARQUEST SCIENTIFIC, INC. IS PROHIBITED.

DO NOT SCALE DRAWING

SIZE	DWG. NO.	REV.
A	LG-DN-R Parts	
SCALE: 1:8	WEIGHT:	SHEET 1 OF 1